

Product Overview

Pfu DNA polymerase catalyzes DNA dependent polymerization from 5'→3', in presence of Mg²⁺ and also exhibits 3'→5' exonuclease activity for proof reading, which enables the polymerase to correct base insertion errors.

Catalog Details

R1201	250 U	Visit us for
R1202	1000 U	more variants
R1203	2000 U	www.dxbidt.com

Storage

- 20 °C

PCR Protocol for templates ≤ 5KB

The samples must be prepared in sterile DNase Free micro-centrifuge tubes, with the following composition placed on ice

Components	50μL reaction	Final Concentration
Template DNA	1 μL	10ng - 250ng
Forward Primer (10μM)	1μL	0.1 - 1μM
Reverse Primer (10μM)	1μL	0.1 - 1μM
dNTPs Mix (2.5mM Each)	4 μL	Upto 200 μM each
10X Pfu Reaction Buffer	5 μL	1X
DX/DT Pfu DNA Polymerase	0.5 μL	0.5 - 1.25 U
Nuclease Free Water	Upto 50 μL	

NOTES:

- Addition of all the above reagents should be done using ice boxes to prevent non-specific amplification.
- To collect all liquid at the bottom of the vial, reaction mix can be kept for a quick spin (10 seconds).
- While doing PCR without heat lid, it is recommended to overlay the sample with 1-2 drops of mineral oil.

Step	Temperature	Time	Cycle
Initial denaturation	95 ± 1 °C	2 - 5 minutes	1
Denaturation	95 ± 1 °C	0.5 - 1 minutes	30 - 35
Annealing *	45 - 65 °C	0.5 - 1 minutes	
Extension**	68 - 72 °C	1 minute/kb	
Final Extension	68 - 72 °C	5-10 minutes	1
Hold, if required.	2- 8 °C	variable	1

PCR Protocol for templates > 5KB

The samples must be prepared in sterile DNase Free micro-centrifuge tubes, with the following composition placed on ice

Components	50μL reaction	Final Concentration
Template DNA	1 μL	50ng- 250ng
Forward Primer (10μM)	1μL	0.1 - 1μM
Reverse Primer (10μM)	1μL	0.1 - 1μM
dNTPs Mix (2.5mM Each)	upto 6 μL	Upto 300 μM each
10X Pfu Reaction Buffer	5 μL	1X
DX/DT Pfu DNA Polymerase	0.5 μL	0.5 - 2.5μL
Nuclease Free Water	Upto 50 μL	

NOTES:

- Addition of all the above reagents should be done using ice boxes to prevent non-specific amplification.
- To collect all liquid at the bottom of the vial, reaction mix can be kept for a quick spin (10 seconds).
- While doing PCR without heat lid, it is recommended to overlay the sample with 1-2 drops of mineral oil.

Step	Temperature	Time	Cycle
Initial denaturation	95 ± 1 °C	2 - 5 minutes	1
Denaturation	95 ± 1 °C	≤ 30s	30 - 35
Annealing *	45 - 65 °C	0.5 - 1 minutes	
Extension**	68 - 72 °C	1.25 minute/kb	
Final Extension	68 - 72 °C	10-20 minutes	1
Hold, if required.	2- 8 °C	variable	1

Quality Control Assays

1. **Purity:** SDS Page analysis with Coomassie Blue Staining resulted in ≥ 95% purity.
2. **Performance testing (1):** 2.5U of enzyme was used to amplify 10ng of DNA template (300bp, 1kb, 3kb and 5kb) in 30 PCR cycles resulted in a single band, confirmed by 1% agarose gel electrophoresis with EtBr.
3. **Performance testing (2):** 2.5U of enzyme was used to amplify 50ng of DNA template (8kb) in 30 PCR cycles resulted in a single band, confirmed by 1% agarose gel electrophoresis with EtBr.
4. **Nuclease tests:** No contamination of endo or exonucleases were detected.

Any Technical Help ?

Please write to us at tech@dxbidt.com . Response can be expected within 24Hrs. Our technical team shall be happy to assist you all the time.

Order Related Queries

Email : info@dxbidt.com | Ph: +91-7349708807

Website: <https://dxbidt.com>

Address: #87, Dasanapura, Lakshmipura Post, Bangalore - 560073