

## Product Overview

Pfu DNA polymerase catalyzes DNA dependent polymerization from 5'→3', in presence of Mg<sup>2+</sup> and also exhibits 3'→5' exonuclease activity for proof reading, which enables the polymerase to correct base insertion errors.

## Catalog Details

R1201	500 U
R1202	2000 U

Visit us for  
more variants  
[www.dxbidt.com](http://www.dxbidt.com)

## Storage

- 20 °C

## PCR Protocol for templates ≤ 5KB

The samples must be prepared in sterile DNase Free micro-centrifuge tubes, with the following composition placed on ice.

Components	50 $\mu$ L reaction	Final Concentration
Template DNA	1 $\mu$ L	<b>10ng</b> - 250ng
Forward Primer (10 $\mu$ M)	1 $\mu$ L	0.1 - 1 $\mu$ M
Reverse Primer (10 $\mu$ M)	1 $\mu$ L	0.1 - 1 $\mu$ M
dNTPs Mix (2.5mM Each)	4 $\mu$ L	Upto 200 $\mu$ M each
10X Pfu Reaction Buffer	5 $\mu$ L	1X
DX/DT Pfu DNA Polymerase	0.5 $\mu$ L	0.5 - 1.25 U
Nuclease Free Water	Upto 50 $\mu$ L	

### NOTES:

- Addition of all the above reagents should be done using ice boxes to prevent non-specific amplification.
- To collect all liquid at the bottom of the vial, reaction mix can be kept for a quick spin (10 seconds).
- While doing PCR without heat lid, it is recommended to overlay the sample with 1-2 drops of mineral oil.

Step	Temperature	Time	Cycle
Initial denaturation	98°C	30s	1
Denaturation	98°C	10s	25-35 Cycles
Annealing *	50 - 65 °C	20s	
Extension**	72 °C	<b>1</b> minute/kb	
Final Extension	72 °C	5-10 minutes	1
Hold, if required.	2- 8 °C	variable	1

## PCR Protocol for templates > 5KB

The samples must be prepared in sterile DNase Free micro-centrifuge tubes, with the following composition placed on ice

Components	50μL reaction	Final Concentration
Template DNA	1 μL	<b>50ng - 250ng</b>
Forward Primer (10μM)	1μL	0.1 - 1μM
Reverse Primer (10μM)	1μL	0.1 - 1μM
dNTPs Mix (2.5mM Each)	upto 6 μL	<b>Upto 300 μM each</b>
10X Pfu Reaction Buffer	5 μL	1X
DX/DT Pfu DNA Polymerase	0.5 μL	<b>0.5 - 2.5μL</b>
Nuclease Free Water	Upto 50 μL	

### NOTES:

- Addition of all the above reagents should be done using **ice boxes** to prevent non-specific amplification.
- To collect all liquid at the bottom of the vial, reaction mix can be kept for a quick spin (10 seconds).
- While doing PCR without heat lid, it is recommended to overlay the sample with 1-2 drops of mineral oil.

Step	Temperature	Time	Cycle
Initial denaturation	98 °C	30s	1
Denaturation	98 °C	10s	25 - 35
Annealing *	50 - 65 °C	20s	
Extension**	72 °C	1 minute/kb	
Final Extension	72 °C	10-20 minutes	1
Hold, if required.	2- 8 °C	variable	1

## Quality Control Assays

1. **Purity:** SDS Page analysis with Coomassie Blue Staining resulted in ≥ 95% purity.
2. **Performance testing (1):** 2.5U of enzyme was used to amplify 10ng of DNA template (300bp, 1kb, 3kb and 5kb) in 30 PCR cycles resulted in a single band, confirmed by 1% agarose gel electrophoresis with EtBr.
3. **Performance testing (2):** 2.5U of enzyme was used to amplify 50ng of DNA template (8kb) in 30 PCR cycles resulted in a single band, confirmed by 1% agarose gel electrophoresis with EtBr.
4. **Nuclease tests:** No contamination of endo or exonucleases were detected.

## CRITICAL NOTE for Site Directed Mutagenesis

Please use 25 to 30 Cycles first to understand amplification results on Gel. Later keep 12- 14 cycles based on yield. Even if products are not visible on gel at low cycles, please go ahead for transformation. Lower the cycles better is fidelity.

## Order Related Queries

Email : [info@dxbidt.com](mailto:info@dxbidt.com) | Ph: +91-7349708807

Website: <https://dxbidt.com>

Address: #87, Dasanapura, Lakshmipura Post, Bangalore - 560073