

## Product Overview

Phi (Φ) Script™ reverse transcriptase is thermostable reverse transcriptase manufactured by dx/dt. It's an engineered form of Moloney Murine Lukemia Virus reverse transcriptase (M-MuLV-RT). Point mutations are introduced in ΦScript™ reverse transcriptase so as to reduce RNase H Activity and to improve its thermal stability. Good cDNA yields can be obtained when reverse transcribed at 42-55°C.

## Catalog Details

*Phi-Script Φ, 200U/μL*

R6101 10,000 Units

R6102 50,000 Units

## Storage

- 20 °C

## First Strand cDNA Synthesis Protocol

The most important element is - good quality RNA. Make sure you have extracted quality RNA (Total RNA or mRNA) from appropriate cells. The quality of the RNA can be checked using agarose gel electrophoresis to see two distinct rRNA. It is important to have intact RNA specially when longer, full length cDNA is expected. Make sure to have all the tips, plastic wares, other consumables and working area RNase Free.

### Full length cDNA synthesis

### Gene Specific cDNA Synthesis

Components	20μL Reaction	Target	20μL Reaction	Target
RNA Template	variable	1 μg	variable	1 μg
dNTP Mix (10mM Each)	1 μL	0.5 mM Each	1 μL	0.5 mM Each
Random Hexamers (50μM Stock) OR Oligo (50μM Each) OR Both	1 μL	upto 200 ng (hexamers) upto 500 ng (oligo-dT)	0 0	NA NA
Gene Specific Primer(10μM Stock) (Reverse ONLY)	0	NA	1 μL	2 - 20pmoles
Total Volume	X μL		X μL	

## PRE - RT STEP

Mix all the above componens ( 1. RNA, 2. dNTP, 3. hexamer OR oligo(dT) or both ) and heat it up at **70°C for 5 minutes**. This step will remove secondary structures of RNA. After this, keep the vial on ice immediately at least for 1 minute. For **gene specific assays**, pre-RT step is optional, may not be essential.

## RT STEP

Perform the following reaction on ice.

Components	20μL reaction	Target Final Concentration
Pre - RT mixture	X μL	
DX/DT Rnase Inhibitor (40U/μL) (Optional)	1 μL	10-40U
5X $\phi$ Script™ Reaction Buffer	4 μL	1X
DTT (100mM Stock)	1 μL	Upto 2 - 10mM
$\phi$ Script™ Reverse Transcriptase (200U/μL)	1 μL	upto 3 μL
Nuclease Free Water	upto 20μL	

- Incubate the above mixture at 42°C for 30 - 60 minutes. Increase the reaction mixture temperature to 50-55°C for **gene specific assays**. If you are using random hexamers, then first incubate 25°C for 10 minutes prior to 42°C.
- Inactivate RT by heating at 85 °C for 10 minutes.
- Store the cDNA samples at -20°C for further analysis. This cDNA now can be used as PCR template for amplification.

## PCR STEP

We recommend to use DX/DT LEO PRIME Mastermix ([R8220](#)) for below PCR assays (End Point PCR ) or Polaris Mastermix ( R2220) for qPCR applications


Components	20μL reaction	50μL reaction	Final Concentration
Template DNA/cDNA	1 μL	2 μL	1ng - 250 ng
Forward Primer (10μM)	0.8 μL	1 μL	0.1 - 1μM
Reverse Primer ( 10μM)	0.8 μL	1 μL	0.1 - 1μM
Extremo MasterMix (2X)	10 μL	25 μL	1X
Nuclease Free Water	Upto 20 μL	Upto 50 μL	

### PCR Conditions

Step	Temperature	Time	Cycle
Initial denaturation	95 °C	2 minutes	1
Denaturation	95 °C	20-30s	25 - 45 Cycles
Annealing (Tm of your primers)	°C	20-30s	
Extension	72 °C	1 minute/kb	
Final Extension	72 °C	2 - 5 minutes	
Hold, if required.	2 - 8 °C	variable	

Any Technical Help ?

Scan the QR Code to read on trouble-shooting guide



CONTACT DETAILS

info@dxbidt.com | www.dxbidt.com | +91-7349708807