

OPTIMUS 100BP DNA LADDER



Product Overview

Optimus™ 100 bp DNA Ladder consists of 12 fragments ranging from 100bp till 1500 bp. The intensity of 500bp and 1000bp are higher compared to other bands and can be used as reference bands. Ladders consists of 100, 200, 300, 400, 500(bright), 600, 700, 800, 900, 1000(bright), 1200, 1500bp. Optimus DNA ladder comes with preloaded dye for ready to use applications.

Catalog Details

R4001-1 0.5ml X 1 vial

R4001-5 0.5ml X 5 vials

Gel % ≥ 1.5

5 μ L per lane loading

100 BP Ladder

FRAGMENT (Conc. per 5 μ L load)

